



# Mind-controlled car demo



## Mind-controlled car demo

- Move your mind-controlled car from the starting line to a given parking region
- The shorter time you use, the higher points you get



## Mind-controlled car

- Use power of alpha band in EEG to control cars
  - Alpha waves increase when we close our eyes and relax
- Use the ratio of alpha power to determine whether the small car move or not
  - If the ratio exceeds a threshold, the car move forward





# Requirements

Only 1 laptop in each group needs to do the following things

1. Clone all source codes (From [Github](#))
2. Install softwares (Follow the instruction in readme)
  - Cygnus
  - Openvibe
  - Python (including the environment)
3. Check the functionality of cars
  - Bluetooth Module
  - Motors

# Requirements - Check the functionality of small car

- Turn on the small car
- Connect the small car to your laptop via bluetooth
  - Find the device called **"CECNL YourGroupNumber"** in the bluetooth setting
  - The password is **1234** for all devices
  - Check the port number





## Requirements - Check the functionality of small car

- Activate the python environment
- Execute a python code “test\_car.py” and enter the port number that the laptop used to connect with the car

```
(lab_car) E:\NYCU PhD\Lab2_car>python test_car.py COM20
```

- Use keyboard (“WASD”) to test the functionality of your car




# Pipeline

1. Put on EEG device and connect to dongle to you laptop
2. Cygnus - Connect to EEG device and show raw EEG signal
3. Openvibe acquisition server - Acquire EEG signal
4. Openvibe designer - Implement real-time EEG processing and calculate alpha power ratio
5. Python - Connect your laptop to small car via bluetooth and send fire command to the car based on current alpha power level

## Cygnus - Connect to EEG device


1. Select the correct device ID on your dongle
2. Connect your device to the laptop







# Cygnus - Check EEG signals

1. Check the impedance of the channels you use
2. Adjust the position of each channel and try to lower the impedance
3. Check the raw EEG signal
  - No too many noise and artifacts



# Openvibe acquisition server

1. Click driver button, and choose LabStreamingLayer
2. Click on driver properties button, and **make sure signal stream is from EEG signal and marker stream is None**
3. Click connect button, and check the terminal
4. Click play button
5. Check the delay time at the bottom of the window. If the delay time is not pretty high, it means that we have successfully received the EEG signals from LSL






# Openvibe designer

1. Open “alpha\_wave\_bandpower.xml”
  - A program for the real-time EEG processing
2. The procedure of EEG processing
  - a. Select a certain EEG channel
  - b. Bandpass filtering
  - c. Extract 4 EEG sub-band power (delta, theta, alpha, beta)
  - d. Calculate the alpha power ratio
3. Send the current alpha power ratio to python via LSL


# Openvibe designer - part 1



- Acquisition client: Receive brain wave data from openvibe acquisition server.
- Channel selector: Double-click on "Channel Selector" and change the "channel list" options to select which channel's signal you want to output .
- Display: Display EEG brain wave signal ( double-click and select "rename" to change the name of signal display to "alpha power")




## Openvibe designer - part 3



Configure Simple DSP settings


Equation  $x*x$

Override settings with configuration file



- Every output of temporal filters are input to “time based epoching” in order to add more sample data. Double click on it, set “duration” to 1 and set “epoch internal” to 0.0625
- Then connect output to “simple DSP”, double click on it and set equation to  $X*X$ , in order to obtain the power of the signal
- Next, output of the “simple DSP” is connected to “signal average” to average the signal

## Openvibe designer - part 5



1. Connect alpha power ratio to “signal display”(named as sign change) to inspect real-time alpha power
2. LSL export (Gipsa): output the alpha power ratio to LSL (we will receive it from python later)
3. Sign change detector: it is used to detect the sign of input data.
4. “Button VRPN Server” and “Button VRPN Client” are use to transfer data
5. Display an image: If the alpha wave exceeds the threshold, the fish picture will be displayed, otherwise it will be empty



## Python part

1. Connect the small car to your laptop via bluetooth
2. Pull the current alpha power ratio from LSL
3. Obtain the average alpha power ratio
4. Set a threshold for average alpha power ratio
  - If it exceed the threshold, send a signal to small car and the car will move forward

```
import serial
import queue
import time

thres = 0.0 # set a threshold here
qsize = 30 # max queue size

def main():
    import argparse
    parser = argparse.ArgumentParser(description="CECNL BCI 2023 Car Demo")
    parser.add_argument("port_num", type=str, help="Arduino bluetooth serial port")
    args = parser.parse_args()
```





## Python part

Assume your computer connects out to car through port “COM(your port number)”

> **test\_car.py**: check car functionalities, you can control the car via keyboard (WASD). To execute, put “COM(your port number)” at the end of argument.

> **mind\_controlled\_car.py**: control car using alpha power ratio of received EEG stream. To execute, put “COM(your port number)” at the end of argument.